




Rocket® Double Lumen Aspiration Needles Reproductive Medicine



- **HIGH HARVESTING & FLUSH RATES.** Rocket® Double Lumen Needles have the highest harvest and flush rates of any comparable needle system.
- **FLUSHING ARM STOPCOCK** fitted as standard to all versions to minimise accidental aspiration of flushing media and increase flushing control.
- **SX™ SPIRAL ECHO BEVEL DESIGN** with high resolution echogenic surface preparation to give smoother, easier tissue penetration with extreme ultrasound definition.
- **SHAPED FINGER GRIP** with bevel marker and enhanced surfacing to significantly increase comfort and control
- **TWO NEEDLE BEVEL STYLES:**
 - Rocket® Embryon®: dual angle bevel with 15mm brushed echogenic band
 - Rocket® SX™: back sharpened bevel, 40mm high definition echo spiral
- **OXFORD TUBING SET:** Single limb tube set into clear silicone bung with 45° metal cannula connection (as illustrated)

Rocket® SX™ Spiral Echo Needles are the latest generation of high performance, high ultrasound definition oocyte aspiration sets from Rocket Medical plc.



Description	Dimensions	Order Code
Rocket® Double Lumen Oocyte Aspiration Needles etched echo band, polyurethane, Oxford style tube set with clear silicone stopper for Falcon® tubes, pouch packed, supplied sterile, 10 per carton, for single patient use	16g x 32cm x 70cm	R57610
	16g x 35cm x 70cm	R57611
	16g x 35cm x 90cm	R57611-00-90
Rocket® SX™ DL Oocyte Aspiration Needle SX™ spiral echo band, polyurethane, Oxford style tube set with clear silicone stopper for Falcon® tubes, pouch packed, supplied sterile, 10 per carton, for single patient use	16G x 33cm x 70cm	R57710-SX-70